

medax

PARASET

SET FOR PARACENTESIS



Paracentesis Set for aspiration of ascetic fluid.

- Waste bag with non-return valve and drainage tap.
- Bags in different capacities: 2 Lt, 5 Lt, 8 Lt, 12 Lt.
- Wide range of needles:
 - Trio of 55 mm needles in the diameters 14 G, 16 G, 19 G;
 - Trio of 80 mm needles in the diameters 14 G, 16 G, 19 G;
 - Veress needle;
 - I.V. cannula with atraumatic tip;
- Connection tubes provided with 3-way stopcock.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

MEDaX

PARASET

SET FOR PARACENTESIS**Ordering information:**

NEEDLES	CAPACITY			
	2 lt	5 lt	8 lt	12 lt
Trio of 55 mm needles	PR55002-00	PR55005-00	PR55008-00	PR55012-00
Trio of 80 mm needles	PR80002-00	PR80005-00	PR80008-00	PR80012-00
Veress Needle	PR00002-00	PR00005-00	PR00008-00	PR00012-00
I.V. introducer needle	PR51002-00	PR51005-00	PR51008-00	PR51012-00
PCS per BOX	50	25	25	15

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
 email: customer-care@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
 Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485:2016 standards.
 In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.
 All products undergo intensive clinical testing and are fully EC and FDA approved.