


medax

MED-I

**BONE MARROW ASPIRATION SYSTEM
WITH ERGONOMIC HANDLE**



EXPLANT/TRANSPLANT VERSION

Bone Marrow aspiration system with ergonomic handle “XI”, also available in the explant/transplant version “EI”.

- Ergonomic handle for a comfortable procedure and a firm grip.
- Triple sharpened tip cannula to allow a non traumatic, fast and easy penetration into the bone marrow cavity.
- Adjustable depth setter.
- Luer Lock connection for syringe.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

BONE MARROW ASPIRATION SYSTEM WITH ERGONOMIC HANDLE

Ordering information:

GAUGE	COLOUR	LENGTH (MM)					
		14 030	14 040	14 050	14 070	14 090	14 110
14	● Green	14 030	14 040	14 050	14 070	14 090	14 110
15	● Blue	15 030	15 040	15 050	15 070	15 090	15 110
16	○ White	16 030	16 040	16 050	16 070	16 090	16 110
18	● Pink	18 030	18 040	18 050	18 070	18 090	18 110

Different sizes on demand

Ordering information:

MED-I ASPIRATION TIP

X	I	GAUGE		LENGTH		0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

MED-I EXPLANT/TRANSPLANT TIP

E	I	GAUGE		LENGTH		0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
 Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

