

**medax**

**MED-J**

**BONE MARROW BIOPSY AND ASPIRATION SYSTEM**



Bone marrow biopsy and aspiration needle to perform bone marrow biopsy and aspiration in the iliac crest.

- Coloured-coded.
- Luer Lock connection for syringe.
- Remover guide to facilitate the sample expulsion.
- Provided with an extra cap.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

**medax**

# MED-J

## BONE MARROW BIOPSY AND ASPIRATION SYSTEM

GAUGE	COLOUR	LENGTH (MM)	
8	● Dark Green	08 100	
11	○ White	11 100	11 150
13	● Light blue	13 050	13 090

Different sizes on demand

### Ordering information:

MED-J



Single box: 10 PCS.

#### Medax Srl Unipersonale

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
email: [customer-care@medax.it](mailto:customer-care@medax.it) • PEC: [medax@legalmail.it](mailto:medax@legalmail.it) • [www.medax.it](http://www.medax.it)

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2016 standards. In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates. All products undergo intensive clinical testing and are fully EC and FDA approved.