

medax

MED-L

**BONE MARROW BIOPSY AND ASPIRATION SYSTEM
WITH SPECIAL EXTRACTION CANNULA**



EXTRACTION CANNULA
"P" VERSION

EXTRACTION CANNULA
"S" VERSION

Bone marrow biopsy and aspiration system provided with special safe and atraumatic extraction cannula.

- Coloured-code ergonomic handle.
- Additional extraction cannula to get the sample in a completely atraumatic and safe way.
- "S" version optimal for multiple procedures.
- Luer Lock connection for syringe.
- Remover guide provided with indicators to facilitate the sample expulsion and to check its length during the procedure.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

MED-L

BONE MARROW BIOPSY AND ASPIRATION SYSTEM WITH SPECIAL EXTRACTION CANNULA

GAUGE	COLOUR	LENGTH (MM)		
8	● Dark Green	08 100	08 150	
11	○ White	11 100	11 150	
13	● Light blue	13 070	13 100	13 150

Med-L "P" version 13G only available on special request
Different sizes on demand

Ordering information:

MED-L "S" VERSION DIAMOND TIP

X	L	GAUGE	LENGTH	0	S

MED-L "P" VERSION FISH MOUTH TIP

X	L	GAUGE	LENGTH	0	P

Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

