

medax

MED-B

BONE MARROW BIOPSY SYSTEM



EXPLANT/TRANSPLANT VERSION

Bone marrow biopsy needle to perform bone marrow biopsy and bone marrow explant/transplant “EM” in the iliac crest.

- Coloured-code ergonomic handle.
- Luer Lock connection for syringe.
- Remover guide provided with indicators to facilitate the sample expulsion and to check its length during the procedure.
- Provided with a special cap for a safe removal of the sample.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

medax

MED-B

BONE MARROW BIOPSY SYSTEM

GAUGE	COLOUR	LENGTH (MM)		
8	● Dark Green	08 100	08 150	
11	○ White	11 100	11 150	
13	● Light blue	13 070	13 100	13 150

Different sizes on demand

Ordering information:

MED-B DIAMOND TIP



MED-B EXPLANT/TRANSPLANT VERSION



Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485:2016 standards.
In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.
All products undergo intensive clinical testing and are fully EC and FDA approved.