

medax

MEDWIRE

REPOSITIONABLE BREAST LOCALIZATION NEEDLE
FOR NON PALPABLE MAMMARY LESIONS



SCREW VERSION



EASY RELEASE VERSION



O-TYPE HOOK WIRE

Repositionable Breast localization needle for non palpable mammary lesions.

- Round (O-type) hook wire.
- Extremely sharp atraumatic tip.
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and sliding stopper.
- Depth indicator on the wire.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

medax

MEDWIRE

REPOSITIONABLE BREAST LOCALIZATION NEEDLE
FOR NON PALPABLE MAMMARY LESIONS

GAUGE Ø	WIRE (MM)	LENGTH (MM)		
20	0,40	20 060	20 090	20 120

Different sizes on demand

Ordering information:

MEDWIRE SCREW VERSION

M	W	GAUGE	LENGTH	1	0
---	---	-------	--------	---	---

MEDWIRE EASY RELEASE VERSION

M	W	GAUGE	LENGTH	2	0
---	---	-------	--------	---	---

Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.