

medax

SONOMED

FULL-CORE BIOPSY SYSTEM
FOR HISTOLOGICAL ASPIRATION BIOPSY



SONOMED, full-core biopsy system for ultrasound guided histological aspiration biopsy.

- Ultra thin sharpened cannula with trocar stylet for a quick and painless percutaneous insertion.
- Syringe of 12ml.
- Compatible with SONOGUN.
- Sterilized by ETO, shelf life 5 years.

medax

SONOGUN

REUSABLE ASPIRATION GUN



SONOGUN, Reusable aspiration gun usable with disposable needle SONOMED.

- Light Weight
- One hand usage.
- Compatible with 10ml and 12ml syringe.

medax
medical devices

Your partner in Biopsy and Special Needles

medax

SONOMED

FULL-CORE BIOPSY SYSTEM FOR HISTOLOGICAL ASPIRATION BIOPSY

GAUGE	COLOUR	LENGTH (MM)	
17	● RED	17 100	17 160
18	● Pink	18 100	18 160
20	● Yellow	20 100	20 160
21	● Light green	21 100	21 160

Different sizes on demand

Ordering information:

SONOMED



Single box: 10 PCS.

SONOGUN



Single box: 1 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485:2016 standards.
In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.
All products undergo intensive clinical testing and are fully EC and FDA approved.