

medax

COX

STANDARD COAXIAL INTRODUCER NEEDLE



Standard coaxial introducer needle used for single access in multiple biopsies.

- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Standard Trocar tip.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

STANDARD COAXIAL INTRODUCER NEEDLE

| BIOPSY NEEDLE GAUGE | COX COLOUR | COX LENGTH (MM) | | | | | | | |
|---------------------------|---------------|-----------------|--------|--------|--------|--------|--------|--------|--------|
| | | 80 | 100 | 120 | 160 | 180 | 200 | 250 | 300 |
| 12 | ● Light green | 11 030 | 11 050 | 11 070 | 11 110 | 11 130 | 11 150 | 11 200 | 11 250 |
| 14 | ● Light blue | 13 030 | 13 050 | 13 070 | 13 110 | 13 130 | 13 150 | 13 200 | 13 250 |
| 16 | ● Blue | 15 030 | 15 050 | 15 070 | 15 110 | 15 130 | 15 150 | 15 200 | 15 250 |
| 18 | ● Red | 17 030 | 17 050 | 17 070 | 17 110 | 17 130 | 17 150 | 17 200 | 17 250 |
| 20 | ● Light brown | 19 030 | 19 050 | 19 070 | 19 110 | 19 130 | 19 150 | 19 200 | 19 250 |

To find the desired code for the coaxial needle, cross the measures of the biopsy needle displayed in the black boxes (if you need a coaxial for a 16G100mm needle, cross the reference to find the correct coaxial code, in this case 15050).
Different sizes on demand

Ordering information:

COX

| | | | | | | | | |
|---|---|-------|--|--------|--|--|---|---|
| C | X | GAUGE | | LENGTH | | | 0 | 0 |
| | | | | | | | | |

Single box: 20 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

