

**medax**

# VERTEBROX

VERTEBROPLASTY NEEDLE



TROCAR TIP

VERTEBROX needle is used in Vertebroplasty procedures to inject cement into the vertebral body.

- Available in metal and plastic handle.
- Ergonomic and radiotransparent handle.
- Centimeter marking on the cannula.
- Low resistance during cement application with smooth internal cannula wall.
- Excellent guiding properties with electropolished cannula surface.
- Available in beveled triple edges sharp tip and trocar tip.
- Luer Lock connection for standard syringes.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

**MEDAX**

# VERTEBROX

## VERTEBROPLASTY NEEDLE

GAUGE	COLOUR	LENGTH (MM)		
		11 100	11 120	11 150
11	○ White	11 100	11 120	11 150
13	● Light blue	13 100	13 120	13 150

Different sizes on demand

### Ordering information:

VERTEBROX WITH METAL HANDLE, BEVELED TIP



VERTEBROX WITH METAL HANDLE, TROCAR TIP



VERTEBROX WITH PLASTIC HANDLE, BEVELED TIP



VERTEBROX WITH PLASTIC HANDLE, TROCAR TIP



Single box: 10 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
 Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2016 standards. In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates. All products undergo intensive clinical testing and are fully EC and FDA approved.