

medax

MED-IO

INTRAOSSSEOUS INFUSION NEEDLE



Intraosseous Infusion needle used during emergency conditions, permitting infusion of drugs and fluids.

- Special handle design for a comfortable procedure and a firm grip.
- Trocar tip to allow a non traumatic, fast and easy penetration.
- Luer Lock connection for syringe.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices





Your partner in Biopsy and Special Needles

medax

MED-IO

INTRASOSSEOUS INFUSION NEEDLE

Ordering information:

| GAUGE | COLOUR | LENGTH (MM) | | |
|-------|---|-------------|--------|--------|
| 14 |  Green | 14 025 | 14 030 | 14 040 |
| 15 |  Blue | 15 025 | 15 030 | 15 040 |
| 16 |  White | 16 025 | 16 030 | 16 040 |
| 18 |  Pink | 18 025 | 18 030 | 18 040 |

Different sizes on demand

Ordering information:

MED-IO



Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customer care@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485:2016 standards.
In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.
All products undergo intensive clinical testing and are fully EC and FDA approved.