

**medax**

# MEDONE ULTRA

PROGRAMMABLE AUTOMATIC BIOPSY SYSTEM



Programmable automatic disposable biopsy system.

- One hand usage.
- 25 mm penetration depth.
- 20 mm sample notch.
- Delay firing mode.
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings.
- Can be used with CT introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years

**medax**  
medical devices

Your partner in Biopsy and Special Needles

**MEDaX**

# MEDONE ULTRA

PROGRAMMABLE AUTOMATIC BIOPSY SYSTEM

GAUGE	COLOUR	LENGTH (MM)				
14	● Green	14 060	14 100	14 120	14 160	14 200
16	○ White	16 060	16 100	16 120	16 160	16 200
18	● Pink	18 060	18 100	18 120	18 160	18 200 18 250
20	● Yellow	20 060	20 100	20 120	20 160	20 200

Different sizes on demand

## Ordering information:

MEDONE ULTRA

M	U	GAUGE	LENGTH	0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

INTRODUCER COAXIAL NEEDLE FOR MEDONE ULTRA

C	T	GAUGE	LENGTH	0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

MEDONE ULTRA WITH COAXIAL NEEDLE IN THE SAME PACKAGE

M	U	GAUGE	LENGTH	C	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

Single box: 10 PCS.

### Medax Srl Unipersonale

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485:2016 standards.  
In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.  
All products undergo intensive clinical testing and are fully EC and FDA approved.