

**medax**

# MEDURO KIT

AUTOMATIC DISPOSABLE BIOPSY SYSTEM  
WITH ANESTHETIC NEEDLE



Automatic disposable spring loaded biopsy system for prostate biopsy.

- One hand usage.
- 22 mm penetration depth.
- Double firing trigger in proximal and distal positions.
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings.

Anesthetic needle.

- 20 Gauge, 200 mm , 250mm & 300mm length.
- Transparent Luer Lock connector.
- Sharpened cannula for a safe and painless insertion and penetration.
- Centimeter markings
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

**medax**

# MEDURO KIT

AUTOMATIC DISPOSABLE BIOPSY SYSTEM WITH ANESTHETIC NEEDLE

BIOPSY NEEDLE			ANESTHETIC NEEDLE		ORDERING SIZE
COLOUR	GAUGE	LENGTH (MM)	GAUGE	LENGTH (MM)	
●	18	200	20	200	18 200
●	18	250	20	250	18 250
●	18	300	20	300	18 300

Different sizes on demand

## Ordering information:

MEDURO KIT

Single box: 10 PCS.



### Medax Srl Unipersonale

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485:2016 standards.  
In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.  
All products undergo intensive clinical testing and are fully EC and FDA approved.